



# A universal, high-performance EKG Signal Processing Engine to reduce clinical burden

Dr Austin Gibbs  
Matthew Fitzpatrick  
Mark Lilburn  
Holly Easlea  
Jonathan Francey

Rebecca Funston  
Jordan Diven  
Stacey Murray  
Oliver GJ Mitchell  
Adrian Condon

Dr Andrew RJ Mitchell  
Dr Benjamin Sanchez  
Dr David Steinhaus

B  S E C U R



# Contents

<b>Abstract</b>	<b>4</b>
<hr/>	
<b>Introduction</b>	<b>6</b>
<hr/>	
<b>Methods</b>	<b>8</b>
HeartKey Algorithm Overview	8
Wet Electrode EKG Data Overview	9
Dry Electrode EKG Data Overview	9
EKG Data Annotation	9
Data Processing	10
Performance Metrics	10
Performance Metrics: Equations	11
<hr/>	
<b>Results &amp; Discussion</b>	<b>12</b>
HeartKey Performance on Wet Electrode Data	12
HeartKey Performance on Motion Based Dry Electrode Data	16
Conclusion	17
<hr/>	
<b>References</b>	<b>18</b>
<hr/>	

# Abstract

Electrocardiogram (EKG) signal conditioning is a vital component in the EKG signal processing chain that ensures effective noise removal and accurate feature extraction.

Signal conditioning aids clinical diagnosis by improving the quality of vital sign metrics, as well as enhancing the performance when deployed on Holters, patches, and wearables.

This study evaluates the performance of the FDA 510(k) cleared HeartKey® Signal Conditioning and QRS peak detection algorithms on a range of annotated public and proprietary EKG databases.<sup>1</sup> To confirm the universal performance of B-Secur's HeartKey Signal Conditioning technology, the datasets were specifically chosen to ensure a challenging variation of real-world representative EKG data, including data from healthy and unhealthy patients, wet and dry electrode types, various lead configurations, hardware sources, and stationary/ ambulatory recordings from clinical and non-clinical settings.

A total of 751 raw EKG files from a broad range of use cases were individually processed through the HeartKey Signal Conditioning algorithm. The algorithm includes several advanced filtering steps to enable significant noise removal and accurate identification of the QRS complex. QRS detection statistics were generated against the annotated EKG files.

HeartKey displayed robust performance across 14 EKG databases (7 public, 7 proprietary), covering a range of healthy and unhealthy patient data, collected with both wet and dry electrodes over multiple lead setups. Over the NSR, MIT-BIH, AHA, and MIT-AF public databases, average QRS sensitivity (Se) and positive predictive values (PPV) of 98.90% and 99.08% were achieved. A similar adaptable performance was observed on noisy

## 14

EKG DATABASES  
(7 PUBLIC, 7 PROPRIETARY)

## 751

RAW EKG FILES  
PROCESSED

## 3,135,366

ANNOTATED BEATS  
ANALYZED

EKG data, with a QRS Se of 93.26% and PPV of 90.53% attained on the challenging NST database. Crucially, HeartKey's performance also translated to the dry electrode space, with an average of 99.22% QRS Se and 99.00% PPV observed on eight dry electrode datasets collected on various hardware sources, including two challenging motion-based collection protocols.

HeartKey demonstrates robust signal conditioning and QRS detection performance across the broad range of tested EKG signals. It should be emphasized that in no way have the algorithms been altered or trained to optimize performance on a given database, meaning that HeartKey is a potential universal solution capable of maintaining a high level of performance across a broad range of clinical and everyday use cases.

# Introduction

Cardiovascular Disease (CVD), an umbrella term encompassing an array of disorders affecting the heart and blood vessels, is the leading cause of death worldwide and a significant burden on global healthcare.<sup>2</sup>


**Early detection and monitoring of CVDs is crucial as it allows the identified conditions to be treated and appropriate medical precautions to be established.**

Due to the wealth of physiological information derived from the heart's electrical signal, electrocardiography is among the most effective diagnostic tools available to aid clinicians in the fight against CVD. Although once restricted to clinical settings, integrating EKG functionality into portable devices allows healthcare professionals to continuously monitor cardiac function remotely over extended periods. The ambulatory approach is compelling and is becoming increasingly valuable in diagnosing and managing cardiac arrhythmias, including atrial fibrillation (AF), which manifest infrequently and inconsistently.<sup>3</sup>

Being able to accurately extract the relevant physiological information from patients amidst the background noise of a non-clinical, unstable environment is vital to ensure no further increases in burden to the clinical pathway.

Electrocardiogram (EKG) signals are characterized by five key features (P, Q, R, S and T waves) pertaining to the direction of electrical signal propagation through the heart at various stages of the cardiac cycle (Figure 1). Variations in characteristic waveform morphology, orientation and frequency can indicate various cardiac conditions, such as arrhythmias and ischemic heart disease, among others. Computer-aided EKG algorithms that process, interpret and autonomously diagnose cardiac abnormalities have emerged as powerful tools to support manual diagnosis by specialists.

The QRS complex, which represents ventricular depolarization, is the most prominent waveform in the EKG and the easiest for algorithms to detect due to its high amplitude.<sup>4</sup> Accurate and reliable algorithmic detection of the QRS complex is crucial. It serves as the basis from which: a) other characteristic waveforms (P & T waves) in the EKG can be identified, and b) critical diagnostic parameters can be derived. The efficient and accurate extraction of the latter is essential as such information acts as the foundation from which more complex algorithms can be constructed.



**Figure 1.** Example EKG waveform highlighting key features

For instance, beat to beat (R-R) intervals are obtained by measuring the time between correctly detected QRS signals and can be used to calculate heart rate (HR), heart rate variability (HRV) and act as an input for Arrhythmia Detection algorithms.

Although clinical recording protocols are standardized, QRS signal morphology can vary significantly from patient to patient.<sup>5,6</sup> To maximize compatibility with QRS detection algorithms, minimizing noise contamination on the EKG signal is essential as it allows the QRS complex to be readily distinguished. Noise contamination can arise from various sources, including 50/60 Hz power line interference, the electrode-skin interface, muscle activity, and general motion artefact noise induced by patient movement. In the context of automatic EKG detection algorithms, noise artefacts are especially problematic as they can trigger false-positive events that obscure valid EKG metrics.

As automated detection algorithms become more common, there is a clear need to input high-quality data to ensure they function to a high performance level. This need is exacerbated in ambulatory monitoring applications, as the levels of noise artefacts produced during daily activities are significantly greater than in a hospital setting.<sup>7</sup> Effective signal conditioning algorithms must be carefully designed to ensure that; a) noise artefacts are not falsely classified as QRS complexes, and b) true QRS complexes are not removed alongside noise during filtering steps. Noise artefacts can obstruct the distinction of true QRS complexes, potentially leading to the missed detection of an important pathological event that can delay or prevent the diagnosis of a cardiac abnormality. Therefore, an effective signal conditioning step must follow the acquisition of raw EKG data to remove excess noise and output EKG signals from which the QRS can be correctly identified, and critical diagnostic parameters obtained.

Despite the plethora of QRS detection algorithms introduced over the past few decades, there still lacks a universal algorithm capable of operating with high accuracy across the wide range of clinically relevant use cases. In this study, we introduce the HeartKey<sup>®</sup> Signal Conditioning and QRS Detection algorithms and evaluate their performance on a total of 14 EKG databases, chosen to represent the inevitable variability in signal quality of real-world EKG data.

Across the broad range of use cases, 3,135,366 annotated beats were analyzed in total. HeartKey demonstrated highly adaptable QRS detection accuracy and positive predictivity in all cases. Due to the low memory footprint and processing requirements of the algorithm, it has the potential to be employed in a host of EKG monitoring applications, both inside and outside of a clinical environment.

# Methods

## HeartKey Algorithm Overview

HeartKey QRS Detection and Heart Rate algorithms are FDA 510(k) cleared as a Class II medical device. Operationally, the algorithms consist of several distinct stages, detailed in Figure 2, designed to produce reliable and robust performance from raw EKG data across a range of EKG lead configurations and hardware sources. The input signal is initially processed through a signal filtering step that has been fine-tuned to operate on signals from various electrode materials and hardware sources. This ensures that the variation in signal quality associated with both methods is adequately dealt with and that the information to be extracted from the processed data is accurate, reliable, and therefore valuable at a clinical level.

After filtering, a clean signal is fed through the QRS detection algorithm. Upon successful identification of the QRS locations, an R-R interval series can be calculated. This feeds into, amongst others, the HR algorithm. Throughout the process, the HeartKey algorithms employ a variety of methods to assess each calculated metric for validity; this ensures robustness and accuracy on even the noisiest of signals. The HeartKey algorithms will also assess and provide an indicator of the EKG signal quality, only outputting a HR value if the signal is deemed to be of sufficient quality.




Figure 2. The flow of EKG data through HeartKey algorithm(s).



## Wet Electrode EKG Data Overview

Five public PhysioNet wet electrode databases were chosen to evaluate the performance of HeartKey QRS Detection and Signal Conditioning algorithms: the NSR database, the MIT-BIH database, the MIT-AF database, the AHA database, and the NST database.<sup>8</sup> The respective databases include 179 raw EKG files, 3,086,647 annotated beats, with a variety of healthy and unhealthy patients and EKG morphologies on both clean and noisy signals. The databases were tested

in line with the AAMI/ANSI EC57 standard<sup>9</sup> – a medical standard that pertains to a “protocol for a reproducible test with clinical requirements and emphasizes the record-by-record presentation of results; that reflect an algorithm’s ability to detect events of clinical significance.” Proprietary wet-electrode data was also collected on an industry gold standard ambulatory device (Bittium Faros 180)<sup>10</sup> using a wet electrode lead II configuration.

## Dry Electrode EKG Data Overview

To demonstrate the applicability of the HeartKey Signal Conditioning and QRS Detection algorithms for integration into the ever-increasing range of EKG-functionalised dry electrode hardware, EKG data originating from a variety of devices and challenging collection protocols, including walking, and running, were used. A database of proprietary dry electrode data was designed to validate QRS detection performance on challenging EKG data. This database represents EKG data collected

from real world use, where EKGs are performed in various settings, through various methods, by various operators, which ultimately results in a significant variation of signal quality. Challenging signals in the database include those with a significant degree of high frequency noise, motion artefacts, low QRS amplitude, irregular rhythms, and variable beat morphologies. Performance on this dry electrode EKG data was evaluated against manual peak annotations.

## EKG Data Annotation

With the exception of MIT-AF database, PhysioNet databases have been independently annotated by cardiologists and the performance of HeartKey was generated against these annotations. Beat annotations for the MIT-AF DB were generated by a minimum of two separate annotators using a computer-based annotation tool. This was followed by a group review of any outstanding annotations, during which highlighted discrepancies were resolved. Proprietary dry electrode databases

were manually annotated by board-certified cardiologists. Databases were annotated individually, followed by a similar group review to ensure agreement on annotations for challenging signals. As manual annotation is the gold standard for EKG performance analysis, these annotations were used as a criterion by which the HeartKey Signal Conditioning algorithm QRS detection performance was compared.

## Data Processing

Each EKG file was individually processed through the HeartKey Signal Conditioning algorithm. QRS Detection results were generated using PhysioNet WFDB programs bxb and sumstats. As standard QRS accuracy measurements employ a wide error window (+/- 150 ms), the precise location of detection within the QRS complex is not important, only that this location remains consistent from beat to beat. Looking at this in isolation could mask variation in where the algorithm picks up the beat. It is therefore beneficial to include HR

accuracy measurements. Many measurements exist for HR, and none are universally accepted. As recommended in Section 4.3.3.1 of ANSI EC57, HR statistic reference annotation files were created for each record, calculating the HR from the reference beat annotations with the same method used in the device. The comparison will generate the Root Mean Square Heart Rate Error to measure the error between the reference and test annotations. Heart Rate Error statistics are generated using WFDB programs mxm and sumstats.

## Performance Metrics

There are four outcomes in which the detector is presented with an input that is either an event or a non-event:

OUTCOME	DESCRIPTION
True Positive (TP)	An event detected correctly.*
False-Negative (FN)	A missed event.*
False-Positive (FP)	A non-event detected as an event.
True Negative (TN)	A non-event correctly rejected.

\*A correctly detected event is defined as a QRS detection location within 150 ms of the QRS annotation, as stated within ANSI (AAMI EC67-2012) standards. If QRS detection is outside the 150 ms window, the beat is missed and classified as a false negative.

### Performance Metrics: Equations

The most common detector performance measures are sensitivity (Se) and positive predictive value (PPV), as detailed below. Sensitivity relates to the ability of the algorithm to identify true events correctly and is calculated using the following equation:

$$Se (\%) = \left( \frac{TP}{TP + FN} \right) \times 100$$

PPV relates to the algorithm's ability to avoid incorrectly detecting false events and is calculated using the following equation:

$$PPV (\%) = \left( \frac{TP}{TP + FP} \right) \times 100$$

Root Mean Square HR Error is calculated as a percentage figure that allows comparison of the HeartKey HR information relative to the reference value:


$$RMSE = \sqrt{\frac{\sum_{i=1}^N (Actual\ HR - Annotated\ HR)^2}{N}}$$

# Results & Discussion

## HeartKey Performance on Wet Electrode Data

Numerous algorithms have reported excellent performance statistics (>99% QRS Se & PPV) on databases for which they have undergone a learning period,<sup>11-13</sup> or when designed to give optimal performance on a given database.<sup>14-16</sup> Without cumbersome training periods, relatively few algorithms have been reported to retain a high performance across multiple databases with different morphologies, cardiac conditions, and signal qualities.<sup>17-19</sup> HeartKey Signal Conditioning and QRS detection algorithms were initially evaluated on five PhysioNet databases: NSR DB, MIT-BIH DB, AHA DB, MIT-AF DB and NST DB.

The majority of EKG data in these databases was collected using a wet electrode modified limb lead II setup. Each database possesses inherent challenges for the HeartKey algorithm to overcome, outlined in Figure 3. The NSR database contained EKG data with high signal quality and no significant arrhythmias and was included to demonstrate the performance of HeartKey on optimal wet electrode EKG data. HeartKey achieved average QRS detection Se and PPV of 99.84% and 99.40% respectively on this healthy dataset.



DATABASE	RECORDING LENGTH	ANNOTATED BEATS	NATURE OF DATA	CHALLENGE OF DATABASE
NSR DB	Channel I only 18 x 1440 min	1,722,008	Healthy patients with no significant arrhythmias	Various beat types
MIT-BIH DB	48 x 30 min	83,978	Hospital patients attending an arrhythmia clinic	Variety of rhythms, morphology, and signal quality. Including complex ventricular, junctional, and supraventricular arrhythmias and conduction abnormalities
AHA DB	78 x 30 min	177,317	Mixture of patients with arrhythmias and normal EKGs	A variety of healthy EKG data, cardiac arrhythmias and other conditions. Includes PVCs, ventricular bi- and trigeminy, ventricular couplets, PVCs, VT, and VF
MIT-AF DB	23 x 600 min	1,081,882	Patients diagnosed with atrial fibrillation	Patients diagnosed with atrial fibrillation (mostly paroxysmal)
NST DB	12 x 30 min	21,462	Artificial noise added to records 118 and 119 of the MIT-BIH DB	Challenging EKGs with varying signal-to-noise levels
Bittium Faros 180	20 x 2.5 min	3,256	Proprietary wet electrode data collected on healthy patients	Gold standard ambulatory Holter monitor

Figure 3. Overview of data from wet electrode databases.

MIT-BIH is a benchmark database of ambulatory EKG recordings containing various arrhythmias and cardiac abnormalities. It is by far the most frequently used database to validate the performance of signal conditioning algorithms in the literature.<sup>20</sup> For this study, Channel I of each recording was analyzed. The 5 minute training period at the beginning of each record was also excluded from analysis as HeartKey does not require a learning phase. The AHA database is another popular public EKG database that contains a range of rhythms, including NSR, alongside numerous less common arrhythmias. Over the two databases, virtually all significant, clinically relevant arrhythmias are covered, ranging from mild conditions, such as tachycardia, to life threatening heart rhythms like ventricular fibrillation.

This gives an excellent indication of HeartKey’s performance on the wide array of real-world clinical EKG data. On both databases, HeartKey achieved >99.60% Average QRS PPV. Average QRS Se performance values of 98.96% and 97.43% were respectively obtained for MIT-BIH and AHA databases (Table 1). When considering the clinical impact 97% QRS Sensitivity would have in the ‘worst case’ scenario (at a maximum HR of 200 bpm), this equates to approximately 2-3 missed or extra beats within a 30s recording.

In relation to arrhythmia detection applications, this would not be deemed to be clinically significant.<sup>21</sup> RMS HR error values of 1.09% and 1.93% were achieved respectively for the MIT-BIH and AHA databases.

Although MIT-BIH and AHA databases do contain ambulatory EKG records, these datasets are among numerous other recording protocols. This means that the average QRS detection statistics are not truly representative of HeartKey performance on ambulatory data, which can be considerably more challenging due to inherent noise contamination through patient movement.

To demonstrate the utility of HeartKey on ambulatory wet electrode EKG data, the MIT-AF database, which contains long term, continuous EKG data from patients diagnosed with atrial fibrillation (AF), was next analyzed. HeartKey achieved Average QRS Se of 99.38% and QRS PPV of 97.37% on this ambulatory dataset. To further highlight the high performance of HeartKey on Holter wet electrode data, a proprietary dataset was collected on healthy subjects using a Bittium Faros 180, a gold standard ambulatory recording device. Similar performance statistics were achieved on this dataset with an average QRS Se and PPV of 99.86% and 99.66%.

MEASUREMENTS	NSR DB	MIT-BIH DB	AHA DB	MIT-AF DB	NST DB	BITTIUM FAROS 180
Detected Beats	1,719,427	82,996	173,117	1,075,071	20,000	3,521
False Positives	10,134	111	618	28,315	2,141	12
False Negatives	2581	982	4200	6,811	1,462	5
QRS Average Sensitivity	99.84	98.96	97.43	99.38	93.26	99.86
QRS Average PPV	99.40	99.86	99.60	97.37	90.53	99.66
RMS HR Error (%)	1.09	1.93	5.42	3.21	19.42	0.61

**Table 1. Performance of HeartKey on wet electrode databases.**

The ability of HeartKey to successfully detect the QRS complex amidst various levels of noise was pushed to its limit using the Noise Stress Test (NST) database. In this public dataset, artificial noise is overlaid on two clean EKG signals from the MIT-BIH database (records 118 and 119), to emulate baseline wander, muscle movement artefacts and electrode motion artefacts. HeartKey achieved a QRS Se of 93.26% and PPV of 90.53% on the NST-DB. Performance metrics for the NST DB, detailed in table 1, are understandably lower than the other four PhysioNet databases, however, HeartKey still performed well even with the extreme presence of noise. Individual records from the NST database were analyzed to determine the level of noise at which the performance of the algorithm becomes affected. Average QRS Se and PPV performance values remain above 99% for SNR levels as low as 12 dB. At an SNR of 6 dB (118e06 + 119e06), the average Se and PPV reduce to 98.26% and 95.18% respectively. At 0 dB and -6 dB, the performance deteriorates markedly when noise becomes equal to or greater than the EKG signal.

Wearable devices with EKG functionality are emerging as powerful tools to detect and remotely monitor cardiac abnormalities outside of a clinical environment.<sup>21</sup> Such devices typically measure the heart's electrical signal using a dry electrode single lead EKG setup.

However, data collection on dry electrode EKG wearables is inherently more challenging for two reasons: the positioning of the wearable device at peripheral locations on the body, such as the wrist or hands, can lead to a reduction in signal amplitude, whereas the increased impedance given by dry electrodes leads to enhanced noise interference. Prior to an effective signal processing step, the combined issues produce data in which the relevant EKG waveforms are buried under noise. The lack of redundancy in single lead EKG set-ups further stresses the need for accurate and reliable signal processing algorithms to extract the

maximum amount of diagnostic information from challenging EKG traces.

To demonstrate the broad applicability of HeartKey Signal Conditioning and QRS Detection algorithms, EKG data was collected on a range of dry electrode devices. The collected datasets contain various challenging single-lead EKG signals, including those with high frequency noise, motion artefacts, low QRS amplitude, irregular beats, and irregular rhythms to represent real-world dry electrode data as closely as possible (Figure 4).

HeartKey QRS detection algorithm displayed strong performance on the range of dry electrode use cases, with an average QRS Se and PPV of 99.13% and 99.00% over the six tested databases (Table 2). The highest QRS detection performance was observed with the chest module, with an average QRS Se of 99.95% and PPV of 99.94%. This is unsurprising as the large dry electrode surface area ensures continuous contact with the skin. As the device is securely strapped in place, there is also less electrode-skin contact movement and therefore fewer noise artefacts. The proximity of the device to the heart also allows a high amplitude QRS to be recorded. This hardware is a stark contrast to prototype wristwatches 1 and 2, which feature dry electrodes with smaller surface areas located at a peripheral site on the body (wrist), resulting in a low amplitude EKG that is more prone to noise contamination. Despite these challenges, HeartKey achieved QRS Se (>99.54%) and PPV (>98.29%) on both tested prototype wristwatches.

The steering wheel similarly suffers due to the proximal location of its data collection (hands/fingers) and as subjects are required to grip the wheel in an unsecured manner, there is inevitable noise contamination from muscle contractions and the moving electrode-skin interface. Average QRS Se and PPV values of 99.29% and 99.07% were achieved despite these difficulties. The most challenging dry

electrode datasets in Table 2 are arguably those collected with the handheld EKG device (Entries 1 & 2). The electrode-skin contact site for these use cases is at the extremity of the body (fingertips), resulting in a low EKG amplitude. To add further difficulty, the datasets were also

collected on unhealthy patients and contain various arrhythmias and ectopic beats. HeartKey achieved QRS Se of 97.98% and 98.14% and QRS PPV of 99.30% and 98.65% on handheld EKG device entries 1 & 2 respectively.

### HeartKey Performance on Dry Electrode Data







HARDWARE	LEAD SETUP	SUBJECT TYPE	RECORDING LENGTH	ANNOTATED BEATS	OVERVIEW	AIMS
 Chest Module	MCLI + Lead I	Healthy	12 x 1 min	4,871	Stress-inducing protocol including mental arithmetic and presentation tasks	Performance during stress tests at variable heart rates
 Wristwatch Prototype 1	Lead I	Healthy	15 x 1 min	983	Stationary recording on wearable prototype wristwatches	Performance on data from low EKG amplitude wristwatches
 Wristwatch Prototype 2	Lead I	Healthy	15 x 1 min	822		
 Steering Wheel	Lead I	Healthy	31 x 1 min	1,848	Stationary recording on steering wheel electrodes wristwatches	Performance on lower quality dry electrode collection device
 Handheld EKG Entry 1	Lead I	Unhealthy	55 x 30-60 s	1,549	Stationary recording of patients with tachycardia, bradycardia, atrial fibrillation	Performance on dry electrode data containing arrhythmias
 Handheld EKG Entry 2	Lead I	Unhealthy	214 x 30-60 s	6,739	Stationary recording of patients with ectopic beats	Performance on dry electrode data containing ectopic beats

Figure 4. Overview of data from dry electrode databases.

HARDWARE	DETECTED BEATS	FALSE POSITIVES	FALSE NEGATIVES	QRS AVERAGE SENSITIVITY (%)	QRS AVERAGE PPV (%)	RMS HR ERROR (%)
Chest Module	4,869	3	2	99.95	99.94	0.63
Wristwatch Prototype 1	978	22	5	99.54	98.29	0.98
Wristwatch Prototype 2	821	13	1	99.87	98.63	1.23
Steering Wheel	1,836	16	12	99.29	99.07	1.86
Handheld EKG Entry 1	1,522	7	217	97.98	99.30	2.00
Handheld EKG Entry 2	6,628	115	111	98.14	98.65	2.57

Table 2. Performance of HeartKey on dry electrode databases.

## HeartKey Performance on Motion Based Dry Electrode Data

Diagnosing intermittent cardiac arrhythmias is a challenge. These arrhythmic episodes can occur infrequently and unpredictably, and they generally require prolonged and repeated cardiac monitoring to be successfully detected.<sup>22</sup> Performing this requires the development of reliable and clinically safe ambulatory monitoring methods. To maximize patient compliance and obtain real-world EKG data, the chosen hardware device needs to be discreet, lightweight, and unobtrusive so the patient can continue as close to as possible an uninterrupted daily life routine. However, this requirement needs to be balanced against maximizing data quality to ensure that false positives and negatives are eliminated as much as possible.

Wearable dry electrode EKG setups are ideal for discreetness and allow patients to maintain regular daily routines. For the EKG functionality within the wearable device to be of clinical value, the algorithms must be capable of maintaining a high performance across the spectrum of motion-based scenarios the patient will enact each day — such as walking or exercising — which will invariably alter EKG signal quality.

The performance of HeartKey algorithms was investigated on internally collected EKG data recorded on consumer-grade chest strap electrodes (MCLI) during two motion-based protocols (Table 3). In the first, subjects were instructed to walk on a treadmill (6 km/h) for 4 minutes. In the second protocol, after an initial warm up period, subjects were asked to run on a treadmill at increasing speeds (subject dependent) over 4 minutes. As expected, the raw EKGs are of poor quality, suffer from frequent baseline wander, and contain large amounts of noise artefacts arising from dry electrode-skin contact movement and patient muscle activity. HeartKey signal processing and QRS detection algorithms performed well on these challenging datasets, achieving Average QRS Se and PPV of 99.09% and 98.15% during the walking protocol (Entry 1), and a QRS Se and PPV of 99.86% and 99.87% during the running protocol.

HARDWARE	ELECTRODE SETUP	PROTOCOL OVERVIEW	RECORDING LENGTH	ANNOTATED BEATS	QRS AVERAGE SENSITIVITY (%)	QRS AVERAGE PPV (%)	RMS HR ERROR (%)
Chest Module Entry 1	Dry MCLI + Lead I	Walking on a treadmill at increasing speeds	21 x 4 min	8,858	99.09	98.15	1.56
Chest Module Entry 2	Dry MCLI + Lead I	Running on a treadmill at increasing speeds	189 x 1 min	19,793	99.86	99.87	0.48

**Table 3.** Performance of HeartKey during motion-based EKG recording protocols.



Although the performance was expected to lower during the running protocol due to the increased levels of motion, the opposite was observed. The buildup of sweat at the electrode-skin interface during the running protocol could explain the results, as this would lead to greater conductance, essentially allowing it to act as a wet electrode.

HeartKey's ability to improve the quality of dry electrode EKG data without interrupting the patient's daily routine will facilitate clinical-grade ambulatory monitoring, hugely enhancing the breadth and depth of available diagnostic hardware and improving their performance both in and out of hospital settings.

## Conclusion

The ability of EKG signal conditioning algorithms to achieve relatively high performance on wet electrode databases for which they have been optimized or trained is not uncommon. However, most EKG signal conditioning algorithms lack universality and are incapable of maintaining the same high level of performance over multiple databases and use cases, limiting their application. We have shown that HeartKey is a crucial tool in pursuing a universal approach to EKG signal conditioning, showing accurate and reliable QRS detection performance across a broad range of clinically relevant datasets. It should again be emphasized that HeartKey algorithms require no learning phases, and in no way have been adjusted to perform better on the tested databases.

As with many clinical EKGs, the data is there, but the signal can be of poor quality and hidden under a range of noise. Prior to the development of high-quality signal conditioning, this data was lost or not actionable. Patients and clinicians could go through multiple repeated investigations, extended periods of monitoring or, in some cases, the condition could be overlooked until a more catastrophic cardiac event occurs. Effective signal conditioning with HeartKey allows clinicians to extract the right data, making the crucial intervention without repeated investigations.



# References

1. HeartKey is a UK Registered Trademark of B-Secur Ltd.
2. Cardiovascular diseases (CVDs). World Health Organisation (WHO).  
**[https://www.who.int/news-room/fact-sheets/detail/cardiovascular-diseases-\(cvds\)](https://www.who.int/news-room/fact-sheets/detail/cardiovascular-diseases-(cvds))**. Last accessed at 15:20 on 12/04/2022.
3. Sana, F., Isselbacher, E. M., Singh, J. P., Heist, E. K., Pathik, B., & Armoundas, A. A. Wearable Devices for Ambulatory Cardiac Monitoring: JACC State-of-the-Art Review. *Journal of the American College of Cardiology*, 2020;75(13):1582.  
**DOI: 10.1016/j.jacc.2020.01.046**
4. do Vale Madeiro, J. P., da Silva Monteiro Filho, J. M., & Rodrigues, P. R. F. Techniques for QRS Complex Detection. *Developments and Applications for ECG Signal Processing: Modeling, Segmentation, and Pattern Recognition*. Academic Press, 2019;(4):89–118.  
**DOI: 10.1016/B978-0-12-814035-2.00010-4**
5. Rijnbeek, P. R., van Herpen, G., Bots, M. L., et al. Normal values of the electrocardiogram for ages 16–90 years. *Journal of Electrocardiology*, 2104;47(6), 914–921. **DOI: 10.1016/j.jelectrocard.2014.07.022**
6. Corrado, D., Biffi, A., Basso, C., Pelliccia, A., & Thiene, G. 12-lead ECG in the athlete: physiological versus pathological abnormalities. *British Journal of Sports Medicine*, 2009;43(9), 669–676. **DOI: 10.1136/bjsm.2008.054759**
7. Kumar, A., Komaragiri, R., & Kumar, M. From Pacemaker to Wearable: Techniques for ECG Detection Systems. *Journal of Medical Systems*, 2018;42(2), 34. **DOI: 10.1007/s10916-017-0886-1**
8. Goldberger, A. L., Amaral, L. A. N., Moody, G. B. PhysioBank, PhysioToolkit, and PhysioNet. *Circulation*, 2000;101(23).  
**DOI: 10.1161/01.cir.101.23.e215**
9. CENELEC - EN 60601-2-27. (2014). Medical electrical equipment - Part 2-27: Particular requirements for the basic safety and essential performance of electrocardiographic monitoring equipment.
10. Bittium. (2018, November 8). Faros Mobile 21 CFR 870.2800 approval letter. U.S. Food and Drug Administration.
11. Xue, Q., Hu, Y. H., & Tompkins, W. J. Neural-network-based adaptive matched filtering for QRS detection. *IEEE Transactions on Biomedical Engineering*, 1992;39(4), 317–329 **DOI: 10.1109/10.126604**
12. Xiang, Y., Lin, Z., & Meng, J. Automatic QRS complex detection using two-level convolutional neural network. *BioMedical Engineering Online*, 2018;17(1), 13. **DOI: 10.1186/s12938-018-0441-4**
13. Cai, W., & Hu, D. RS Complex Detection Using Novel Deep Learning Neural Networks. *IEEE Access*, 2008;8, 97082–97089.  
**DOI: 10.1109/ACCESS.2020.2997473**
14. Pan, J., & Tompkins, W. J. A Real-Time QRS Detection Algorithm. *IEEE Transactions on Biomedical Engineering*, 1985;32(3), 230–236. **DOI: 10.1109/TBME.1985.325532**
15. Rahul, J., Sora, M., & Sharma, L. D. Dynamic thresholding based efficient QRS complex detection with low computational overhead. *Biomedical Signal Processing and Control*, 2021;67, 102519. **DOI: 10.1016/j.bspc.2021.102519**
16. Farashi, S. A multiresolution time-dependent entropy method for QRS complex detection. *Biomedical Signal Processing and Control*, 2016;24, 63–71. **DOI: 10.1016/j.bspc.2015.09.008**
17. Kunzmann, U., von Wagner, G., Schöchlin, J., & Bolz, A. Parameter extraction of ECG signals in real-time. *Biomedizinische Technik. Biomedical Engineering*, 2002;47 (2), 875–878. **DOI: 10.1515/bmte.2002.47.s1b.875**
18. Dotsinky, I. A., & Stoyanov, T. v. Ventricular beat detection in single channel electrocardiograms. *BioMedical Engineering Online*, 2004;3(1), 1–9. **DOI: 10.1186/1475-925X-3-3**
19. Kim, J., & Shin, H. Simple and Robust Realtime QRS Detection Algorithm Based on Spatiotemporal Characteristic of the QRS Complex. *PLOS ONE*, 2016;11(3), e0150144. **DOI: 10.1371/journal.pone.0150144**
20. Moody, G. B., & Mark, R. G. The impact of the MIT-BIH Arrhythmia Database. *IEEE Engineering in Medicine and Biology Magazine*, 2001;20(3), 45–50. **DOI: 10.1109/51.932724**
21. Bouzid, Z., Al-Zaiti, S. S., Bond, R., & Sejdic, E. Remote and wearable ECG devices with diagnostic abilities in adults: A state-of-the-science scoping review. *Heart Rhythm*. In press. **DOI: 10.1016/j.hrthm.2022.02.030**
22. Heidt, S. T., Kratz, A., Gonzalez, R., et al. Symptoms In Atrial Fibrillation: A Contemporary Review And Future Directions. *Journal of Atrial Fibrillation*, 2016;9(1), 1422. **DOI: 10.4022/jafib.1422**





Find out more about  
us on our website

[info@b-secur.com](mailto:info@b-secur.com)

[in](#) @b-secur

[@bsecurtech](#)